


Product Datasheet

Mouse Pro-MMP9 ELISA Kit (orb1289395)

Description	Mouse Pro-MMP9 ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	MMP9
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P41245
Sensitivity	156 pg/ml
Entrez	17395
Expiration Date	Please enquire.



Representative
Standard
Curve