


Product Datasheet

Mouse IgG ELISA Kit (orb1289408)

Description	Mouse IgG ELISA Kit
Reactivity	Mouse
Range	1560 pg/ml - 100000 pg/ml
Tested Applications	ELISA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Sensitivity	800 pg/ml
Expiration Date	Please enquire.



Representative
Standard
Curve