


Product Datasheet

Mouse IFN beta ELISA Kit (orb1289412)

Description	Mouse IFN beta ELISA Kit
Reactivity	Mouse
Range	12.5 pg/ml - 800 pg/ml
Tested Applications	ELISA
Target	IFNB1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P01575
Sensitivity	6 pg/ml
Entrez	15977
Expiration Date	Please enquire.



Representative
Standard
Curve