


Product Datasheet

Mouse HSP70 ELISA Kit (orb1289413)

Description	Mouse HSP70 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	HSPA1A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	Q61696
Sensitivity	80 pg/ml
Entrez	193740
Expiration Date	Please enquire.



Representative
Standard
Curve