


Product Datasheet

Mouse CXCL7 ELISA Kit (orb1289419)

Description	Mouse CXCL7 ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	PPBP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9EQI5
Sensitivity	4 pg/ml
Entrez	57349
Expiration Date	Please enquire.



Representative
Standard
Curve