


Product Datasheet

Mouse CXCL12 ELISA Kit (orb1289421)

Description	Mouse CXCL12 ELISA Kit
Reactivity	Mouse
Range	46.9 pg/ml - 3000 pg/ml
Tested Applications	ELISA
Target	CXCL12
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P40224
Sensitivity	25 pg/ml
Entrez	20315
Expiration Date	Please enquire.



Representative
Standard
Curve