


Product Datasheet

Mouse CHL1 ELISA Kit (orb1289423)

Description	Mouse CHL1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	CHL1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P70232
Sensitivity	80 pg/ml
Entrez	12661
Expiration Date	Please enquire.



Representative
Standard
Curve