


Product Datasheet

Mouse CD258 ELISA Kit (orb1289424)

Description	Mouse CD258 ELISA Kit
Reactivity	Mouse
Range	93.7 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	TNFSF14
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	Q9QYH9
Sensitivity	50 pg/ml
Entrez	50930
Expiration Date	Please enquire.



Representative
Standard
Curve