


Product Datasheet

Mouse CCL4 ELISA Kit (orb1289427)

Description	Mouse CCL4 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CCL4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P14097
Sensitivity	7 pg/ml
Entrez	20303
Expiration Date	Please enquire.



Representative
Standard
Curve