


Product Datasheet

Human ULBP2 ELISA Kit (orb1289434)

Description	Human ULBP2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	ULBP2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	Q9BZM5
Sensitivity	30 pg/ml
Entrez	80328
Expiration Date	Please enquire.



Representative
Standard
Curve