


Product Datasheet

Human Ubiquitin ELISA Kit (orb1289435)

Description	Human Ubiquitin ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	UBB
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P0CG47
Sensitivity	30 pg/ml
Entrez	7314
Expiration Date	Please enquire.



Representative
Standard
Curve