


Product Datasheet

Human Peroxiredoxin 2 ELISA Kit (orb1289447)

Description	Human Peroxiredoxin 2 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	PRDX2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P32119
Sensitivity	156 pg/ml
Entrez	7001
Expiration Date	Please enquire.



Representative
Standard
Curve