


Product Datasheet

Human MMP10 ELISA Kit (orb1289452)

Description	Human MMP10 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	MMP10
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P09238
Sensitivity	25 pg/ml
Entrez	4319
Expiration Date	Please enquire.



Representative
Standard
Curve