


Product Datasheet

Human Ki67 ELISA Kit (orb1289456)

Description	Human Ki67 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	MKI67
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P46013
Sensitivity	8 pg/ml
Entrez	4288
Expiration Date	Please enquire.



Representative
Standard
Curve