


Product Datasheet

Human IL-1F5 ELISA Kit (orb1289457)

Description	Human IL-1F5 ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	IL36RN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q9UBH0
Sensitivity	50 pg/ml
Entrez	26525
Expiration Date	Please enquire.



Representative
Standard
Curve