


Product Datasheet

Human HIF1 alpha ELISA Kit (orb1289468)

Description	Human HIF1 alpha ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	HIF1A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q16665
Sensitivity	60 pg/ml
Entrez	3091
Expiration Date	Please enquire.



Representative
Standard
Curve