


Product Datasheet

Human FGFR3 ELISA Kit (orb1289473)

Description	Human FGFR3 ELISA Kit
Reactivity	Human
Range	250 pg/ml - 16000 pg/ml
Tested Applications	ELISA
Target	FGFR3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P22607
Sensitivity	125 pg/ml
Entrez	2261
Expiration Date	Please enquire.



Representative
Standard
Curve