


Product Datasheet

Human FGFBP1 ELISA Kit (orb1289474)

Description	Human FGFBP1 ELISA Kit
Reactivity	Human
Range	23.4 pg/ml - 1500 pg/ml
Tested Applications	ELISA
Target	FGFBP1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q14512
Sensitivity	10 pg/ml
Entrez	9982
Expiration Date	Please enquire.



Representative
Standard
Curve