


Product Datasheet

Human CD170 ELISA Kit (orb1289487)

Description	Human CD170 ELISA Kit
Reactivity	Human
Range	46.9 pg/ml - 3000 pg/ml
Tested Applications	ELISA
Target	SIGLEC5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	O15389
Sensitivity	25 pg/ml
Entrez	8778
Expiration Date	Please enquire.



Representative
Standard
Curve