


Product Datasheet

Human BMP9 ELISA Kit (orb1289489)

Description	Human BMP9 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	GDF2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Bone Tissue Homogenates, Cell Culture Supernates, Serum
Uniprot ID	Q9UK05
Sensitivity	7 pg/ml
Entrez	2658
Expiration Date	Please enquire.



Representative
Standard
Curve