


Product Datasheet

Human FGF4 ELISA Kit (orb1289499)

Description	Human FGF4 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	FGF4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P08620
Sensitivity	10 pg/ml
Entrez	2249
Expiration Date	Please enquire.



Representative
Standard
Curve