


Product Datasheet

ABHD14B antibody (orb128950)

Description	Rabbit Polyclonal to ABHD14B.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96IU4
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...