

Product Datasheet

Sheep IgG protein (orb1289502)

Description	Sheep IgG protein
Tested Applications	IH, IP, WB
Preservatives	Purified Sheep IgG in PBS, pH 7.3, 0.01% sodium azide.
Form/Appearance	Purified Sheep IgG in PBS, pH 7.3, 0.01% sodium azide.
Storage	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Please minimize freeze-thaw cycles.
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Sheep
Expiration Date	6 months from date of receipt.