

Product Datasheet

COVID-19 S-trimer Protein (G614&K439) (orb1289530)

Description	Recombinant COVID-19 S-trimer Protein (G614&K439)
Reactivity	Virus
Tested Applications	ELISA, MS, SDS-PAGE, WB
Preservatives	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Please minimize freeze-thaw cycles.
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Human cells
Uniprot ID	P0DTC2
Entrez	43416, 1
Expiration Date	6 months from date of receipt.