

Product Datasheet

EPHA2/3/4 antibody (orb128954)

Description	Rabbit polyclonal antibody to EPHA2/3/4
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P29317, P29320, P54764
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.

