



Product Datasheet

BAF250B antibody (orb128959)

Description	Rabbit Polyclonal to BAF250B.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NFD5
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:00-1:500, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from L...



Confocal immunofluorescence analysis of ...