

Product Datasheet

Rabbit Anti-Chicken IgG (H&L) antibody (Biotin) (orb1289675)

Description	Rabbit Polyclonal antibody against Rabbit Anti-Chicken IgG conjugated to Biotin
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, IH, WB
Immunogen	Chicken IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	E: 1:5000-1:20000, WB: 1:5000-1:20000, IH: 1:500-1:1000, IC: 1:500-1:1000
Expiration Date	12 months from date of receipt.