

Product Datasheet

Mouse Anti-Human IgG4 antibody (orb1289698)

Description	Mouse Monoclonal antibody to Mouse Anti-Human IgG4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IH, WB
Immunogen	Human IgG4
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	E: 1:5000-1:20000, WB: 1:5000-1:20000, IH: 1:500-1:1000
Expiration Date	12 months from date of receipt.