



Product Datasheet

NEFM antibody (orb1289733)

Description	Rabbit Polyclonal antibody to NEFM
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human of human NEFM. The exact sequence is proprietary.
Target	NEFM
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P07197 , P08553 , P12839
Entrez	4741 , 24588 , 18040
Dilution Range	WB: 1:500-1:1000, IH: 1:50-1:200, IF:IC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NEFM expression...



Immunofluorescent analysis of NEFM stain...