



Product Datasheet

CD107a antibody (orb1289743)


Description	Rabbit Polyclonal antibody to CD107a
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human CD107a. The exact sequence is proprietary.
Target	LAMP1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P11438 , P11279 , P14562
Entrez	25328 , 3916 , 16783
Dilution Range	WB: 1:500-1:1000, IH: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of CD107a expressi...



Immunohistochemical analysis of CD107a
S...



Immunofluorescent analysis of CD107a
sta...