



Product Datasheet

Arginase 1 antibody (orb1289744)

Description	Rabbit Polyclonal antibody to Arginase 1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human Arginase 1. The exact sequence is proprietary.
Target	ARG1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P07824 , P05089 , Q61176
Entrez	11846 , 383 , 29221
Dilution Range	WB: 1:500-1:1000, IH: 1:50-1:100, IF:IC: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of Arginase 1 expr...



Immunofluorescent analysis of Arginase 1...