


Product Datasheet

Histone H1 antibody (orb1289767)

Description	Rabbit Polyclonal antibody to Histone H1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Histone H1. The exact sequence is proprietary.
Target	HIST1H1A; HIST1H1C; HIST1H1D; HIST1H1E
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P43275 , P16403 , Q02539 , P10412 , P15864 , D4A3K5 , P43277 , P43274 , P15865 , P16402
Entrez	291145 , 3007 , 3006 , 3008 , 50708 , 80838 , 50709 , 3024 , 14957 , 201097
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H1 expr...