

Product Datasheet

Histone H4 (MonoMethyl-R3) antibody (orb1289771)

Description

Rabbit Polyclonal antibody to Histone H4

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

KLH-conjugated synthetic Mono-methylated peptide corresponding to residues surrounding R3 of human Histone H4 protein. The exact sequence is proprietary.

Target

HIST1H4A; HIST1H4B; HIST1H4C; HIST1H4D; HIST1H4E; HIST1H4F; HIST1H4H; HIST1H4I; HIST1H4J; HIST1H4K; HIST1H4L; HIST2H4A; HIST2H4B; HIST4H4

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P62806, P62804, P62805

Entrez

8361, 319159, 102557184, 100041230, 319158, 8360, 69386, 502913, 8370, 8368, 8294, 102548682, 319156, 100360950, 121504, 319161, 8359, 8366, 500351, 319157, 295277, 320332, 319155, 8367, 8365, 8364, 291152, 100912290, 102641229, 8363, 554313, 102551184, 326619, 319160, 680097, 97122, 100912564, 64627, 100912418, 8362, 326620


Dilution Range

WB: 1:500-1:1000


Expiration Date

12 months from date of receipt.

kDa A B C



Western blot analysis of Histone H4 (Mon...)



Direct ELISA antibody dose-response curve...