

Product Datasheet

SLC22A6 peptide (orb12898)

Description	SLC22A6 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	SLC22A6
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb11177 (SLC22A6 antibody)
Expiration Date	6 months from date of receipt.