



Product Datasheet

DARPP32 antibody (orb1289866)

Description	Rabbit Polyclonal antibody to DARPP32
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human DARPP32. The exact sequence is proprietary.
Target	PPP1R1B
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9UD71 , Q60829 , Q6J4I0
Entrez	19049 , 84152 , 360616
Dilution Range	WB: 1:500-1:1000, IH: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of DARPP32 express...



Immunohistochemical analysis of DARPP32

...