

Product Datasheet

Histone H4 antibody (orb1289888)

Description

Rabbit Polyclonal antibody to Histone H4

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

IH, WB

Applications
Immunogen

KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Histone H4. The exact sequence is proprietary.

Target

HIST1H4A; HIST1H4B; HIST1H4C; HIST1H4D; HIST1H4E; HIST1H4F; HIST1H4H; HIST1H4I; HIST1H4J; HIST1H4K; HIST1H4L; HIST2H4A; HIST2H4B; HIST4H4

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source


Rabbit

Uniprot ID
P62806, P62804, P62805
Entrez
8361, 319159, 102557184, 100041230, 319158, 8360, 69386, 502913, 8370, 8368, 8294, 102548682, 319156, 100360950, 121504, 319161, 8359, 8366, 500351, 319157, 295277, 320332, 319155, 8367, 8365, 8364, 291152, 100912290, 102641229, 8363, 554313, 102551184, 326619, 319160, 680097, 97122, 100912564, 64627, 100912418, 8362, 326620
Dilution Range


WB: 1:500-1:1000, IH: 1:50-1:200

Expiration Date

12 months from date of receipt.

kDa A B C D


Western blot analysis of Histone H4 expr...



Immunohistochemical analysis of Histone ...