


Product Datasheet

S100A16 antibody (orb128989)

Description	Rabbit Polyclonal to S100A16.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96FQ6
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 3...