

Product Datasheet

Histone H4 (MonoMethyl-K31) antibody (orb1289903)

Description

Rabbit Polyclonal antibody to Histone H4

kDa A B C
70

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Western blot analysis of Histone H4 (Mon...

Immunogen

KLH-conjugated synthetic Mono-methylated peptide corresponding to residues surrounding K31 of human Histone H4 protein. The exact sequence is proprietary.

Target

HIST1H4A; HIST1H4B; HIST1H4C; HIST1H4D; HIST1H4E; HIST1H4F; HIST1H4H; HIST1H4I; HIST1H4J; HIST1H4K; HIST1H4L; HIST2H4A; HIST2H4B; HIST4H4

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P62806, P62804, P62805

Entrez

8361, 319159, 102557184, 100041230, 319158, 8360, 69386, 502913, 8370, 8368, 8294, 102548682, 319156, 100360950, 121504, 319161, 8359, 8366, 500351, 319157, 295277, 320332, 319155, 8367, 8365, 8364, 291152, 100912290, 102641229, 8363, 554313, 102551184, 326619, 319160, 680097, 97122, 100912564, 64627, 100912418, 8362, 326620

Dilution Range

WB: 1:500-1:1000

Expiration Date

12 months from date of receipt.