


Product Datasheet

Syndecan 1 antibody (orb128991)

Description	Rabbit polyclonal antibody to Syndecan 1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P18827
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
extracts
from H...