

Product Datasheet

IgG control Antibody (orb1289952)

Description	Mouse IgG2b Kappa Isotype Control (mAb)
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG control
Preservatives	Powder
Form/Appearance	Liquid
Storage	-20°C to -70°C for 12 months in lyophilized state. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Note	For research use only
Application notes	This antibody is suitable for use as a non-targeting isotype control in various in vitro and in vivo studies. It can also be used as a negative control in various applications such as ELISA, Western blot, immunofluorescence, immunohistochemistry, immunoprecipitation, and flow cytometry. Each laboratory should determine an optimum working titer for use in its particular application.
Isotype	Mouse IgG2b Mouse kappa
Clonality	Recombinant
Purity	95%
Expiration Date	12 months from date of receipt.