


Product Datasheet

Human IL-22 R alpha 1 Protein (orb1289953)

Description	Interleukin-22 receptor subunit alpha-1(IL22RA1) is a component of the receptor for IL20, IL22 and I...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-22 R alpha 1
Preservatives	51 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 50.9 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-22 R alpha 1, Fc Tag (orb1289953) is expressed from human 293 cells (HEK293). It contains AA His 16 - Thr 228 (Accession # Q8N6P7-1).
MW	50.9 kDa
NCBI	NP_067081.2
Expiration Date	6 months from date of receipt.



Loaded
Human IL-
22 R alpha
1 Protein,
Fc...