

Product Datasheet

Biotinylated Human IL-18 R beta / IL-1 R7 Protein (orb1289955)

Description

Biotinylated Human IL-18 R beta,

cells (HEK293). It contains AA Phe...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

IL-18 R beta

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 42.0 kDa. The protein migrates as 55-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Biotinylated Human IL-18 R beta, His, (orb1289955) is expressed from human 293 cells (HEK293). It contains AA Phe 20 - Arg 356 (Accession # O95256-1).

MW

42.0 kDa

Uniprot ID

O95256

Expiration Date

6 months from date of receipt.