

Product Datasheet

Fibroblast activation protein, alpha/FAP Antibody (orb1290013)

Description

Fibroblast activation protein, alpha/FAP Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, WB

Applications

Immunogen

E.coli-derived human Fibroblast activation protein, alpha/FAP recombinant protein (Position: Q65-Q580).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

97 kDa

Uniprot ID


[Q12884](#)

Dilution Range


Western blot, 0.1-0.25 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human

Expiration Date

12 months from date of receipt.



Western blot analysis of Fibroblast acti...



Flow Cytometry analysis of SiHa cells us...