


Product Datasheet

Human eNOS/NOS3 ELISA Kit (orb1290034)

Description	Human eNOS/NOS3 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G2-P1203
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants
Uniprot ID	P29474
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in cell culture supernatants is less than the lowest standard, 78 pg/ml.
Expiration Date	Please enquire.



Human
NOS3
PicoKine
ELISA Kit
standard c...