

Product Datasheet

dsRNA antibody (orb1290052)

Description	Mouse monoclonal antibody to dsRNA
Species/Host	Mouse
Endotoxins	Please contact with the lab for this information.
Conjugation	Unconjugated
Tested Applications	DOT, ELISA, IF
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG2a, κ
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.