

Product Datasheet

ACE2 mouse antibody (orb1290117)

Description	Mouse monoclonal antibody to ACE2 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Human ACE2 protein
Preservatives	PBS, pH 7.4, containing 0.05% proclin300, 50% glycerol.
Form/Appearance	Liquid
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -99 °C for twelve months from the date of receipt.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	WB:1:1000-4000;IHC:1:50-100
Expiration Date	12 months from date of receipt.