

Product Datasheet

SARS-CoV-2(2019-nCoV) Spike glycoprotein antibody (orb1290118)

Description	Mouse monoclonal antibody to SARS-CoV-2(2019-nCoV) Spike glycoprotein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Spike protein fragment 2[Y1047-Q1210]
Preservatives	PBS, pH 7.4, containing 0.05% proclin300, 50% glycerol.
Form/Appearance	Liquid
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -99 °C for twelve months from the date of receipt.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	ELISA:1:4000-1:8000,WB:1:1000-5000
Expiration Date	12 months from date of receipt.