

Product Datasheet

SARS-CoV-2(2019-nCoV) Nucleocapsid antibody (orb1290138)

Description	Human monoclonal antibody to SARS-CoV-2(2019-nCoV) Nucleocapsid
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Preservatives	PBS,pH7.5
Form/Appearance	Liquid
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only
Isotype	Human IgG
Clonality	Monoclonal
Dilution Range	WB:1:1000-2000; ELISA:1:5000- 1:10000
Expiration Date	12 months from date of receipt.