

Product Datasheet

Human ACE2 protein (orb1290161)

Description	Recombinant Human ACE2 protein
Reactivity	Human
Endotoxins	Please contact with the lab for this information
Conjugation	Unconjugated
Preservatives	Supplied as solution form in TBS,pH8.0/ Supplied as lyophilized from TBS,pH8.0
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only
Purity	> 90% as determined by SDS-PAGE
MW	85.2kDa
Expiration Date	6 months from date of receipt.