

Product Datasheet

SARS-CoV-2(2019-nCoV) 3C-like Proteinase Protein (orb1290165)

Description	Recombinant SARS-CoV-2(2019-nCoV) 3C-like Proteinase Protein
Reactivity	Virus
Endotoxins	Please contact with the lab for this information
Conjugation	Unconjugated
Preservatives	Supplied as solution form in PBS, pH 7.5/ Supplied as lyophilized from PBS, pH 7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only
Purity	> 90% as determined by SDS-PAGE
MW	34.95kDa
Expiration Date	6 months from date of receipt.